



The Walking Pants can only be applied as depicted above so that the Loop Straps 1a and 1b are at the back and over the coccyx (tailbone).

Product Features

- The Walking Pants are designed as an aid to assist a patient who has difficulty in walking. The Walking Pants may also be used in rehabilitation setting for gait training.
- Walking Pants are generally designed for use with a ceiling lift with a track system. Ceiling tracks offer more versatility and freedom of movement. Because of height restrictions, the use with a mobile floor lift may not be suitable for taller patients.
- Please note that the Pants are intended for use with patients who are able to actively participate in achieving a standing position and are able to support their weight.

Applying the Walking Pants

The Walking Pants should not be applied with the patient standing. They should be applied with the patient sitting.

Seated Patient On Bed

- Seat the patient to the edge of the bed with one caregiver on either side of the patient.
- Lean the patient to one side and into the caregiver. The caregiver should lean the patient to such an extent so that the patient's buttocks on the opposite side are slightly off the bed.
- The caregiver on the opposite side, should slide the Walking Pants under the lifted buttocks.
- Repeat for the opposite side and with a slightly rocking motion, adjust the Walking Pants so that patient is in the middle and center of the pants.
- Attach the loops and straps as described below.

Seated Patient in Chair

- Follow the same procedures as applying a Pants making sure that the Pant with loop straps 1a and 1b are at the back. Attach the loops and straps as described below

