



TONYA – Child's lift
The special lift for children





TONYA –
motorized sitting/lying frame

TONYA

The special child's lift

Horcher's TONYA child's lift has been developed from Horcher's two most successful original lift products (DIANA and LEXA) and therefore combines functionality, design and exceptional quality. Of highly colorful design, TONYA permits children with movement impairments to be lifted from any position in an ergonomically correct manner.

Other features of the product include:

- › The legs of the lift can be opened to two positions to enable the patient to be lifted from a wheel chair, for example;
- › Infinitely adjustable and jolt-free adjustment of the lifting column;
- › Electric lift drive;
- › IBS – Intelligent Battery System;
- › Jolt-free and safe movement without strain thanks to smooth-running castors;
- › Locking brakes on the wheels;
- › Wide range of slings for a diversity of uses;
- › Lifting capacity of 440 pounds



With the motorized sitting/lying frame the child can be effortlessly moved without any physical strain from a sitting to a lying position, for example.

Your distributor:

**BARRIER
FREE LIFTS**

Barrier Free Lifts™
1620 SW 17th St
Ocala, Florida 34471
USA

Phone: 800-582-8732
Fax: 866-378-3318
www.barrierfreelifts.com
info@barrierfreelifts.com