



***LEXA – Professional Lift***  
***The multi-talent among lifts***





Sitting position ...



Lying position ...



Digital scale as LEXA accessory

## LEXA

### A multi-talent for careful transport of patients

The Original LEXA lift is an absolute multi-talent among lifts. It is a very comfortable belt lifter, which, due to its swivel spreader bar, has many capabilities and applications. For instance, it can easily take up patients sitting, lying, or from the floor and move them easily and safely to a sitting or lying position – so they can be transported gently and comfortably to the toilet, to the wheelchair, or to the bed. The sway-free transport also conveys a feeling of security. The LEXA lift is available with a wide variety of slings; it can be used both with loop slings and with 4-point clip connectors.

- › Spreadable base in three positions to receive the patient from the wheelchair;
- › Electric lift motor;
- › IBS – Intelligent Battery System;
- › Effortless, smooth, safe transport thanks to smooth-running casters;
- › Lockable casters;
- › Versatile sling program for the widest variety of applications;
- › Lifting capacity 440 pounds



LEXA with motorized spreader bar



Dolphin patient positioner



Swivel bar

Your distributor:

**BARRIER  
FREE LIFTS**

Barrier Free Lifts™  
1620 SW 17th St  
Ocala, Florida 34471  
USA

Phone: 800-582-8732  
Fax: 866-378-3318  
[www.barrierfreelifts.com](http://www.barrierfreelifts.com)  
[info@barrierfreelifts.com](mailto:info@barrierfreelifts.com)