



The DIANA patient lift is an ideal aid for transfer to the bath, to the toilet, or wherever the patient would like to go. The casters ensure a smooth transport and can be conveniently block by parking brakes. The lifting rod works reliably, continuously, as well as smooth and granted patients the highest possible security.

The transfer lift DIANA can help to make transfers safer and easier for you and the person that you are moving. The DIANA sling lift has for years been a guarantee for qualified care and satisfied patients in nursing homes and home care. For persons with limited mobility, it provides ergonomically correct lifting from any position.

Features at a glance

Patient Lift DIANA Comfort

- › Spreadable base e.g. for use when lifting from a wheelchair
- › Smooth and swayfree lifting via hand controls
- › Removable 24 V / 2.9 Ah Battery-Pack incl. wall charging device
- › Spreader bar with 4-point loop connector
- › Emergency Off button
- › Manual and electric emergency lowering
- › Safe, effortless transport by 4 light casters 100 mm
- › Lockable brakes
- › Comprehensive range of slings to meet your individual requirements
- › Lifting capacity up to 440 pounds (lbs.)
- › Color: Perl Gold / Anthracite
- › Optional: integrated digital patient scale



Digital Scale (optional)



4-Point Spreader Bar

Base length	1187 mm
Base width closed/open	683/1043 mm
Base height	137 mm
Overall height min./max.	1446 mm / 2013 mm
Lifting range	989 mm
Spreader bar	4-point loop connector
Lifting capacity	up to 440 pounds (lbs.)
Battery capacity	24 V / 2.9 amps