

Patient Lift LEXA with motorized patient positioner

Features at a glance

Professional Patient Lift LEXA with motorized patient positioner

- › Motorized patient positioner with clips suspension adjustable via hand control e.g. to bring a person continuously into a sitting or lying position
- › Spreadable base e.g. for use when lifting from or into a wheelchair
- › Smooth and swayfree lifting via hand controls
- › Safe, effortless transport by 4 light casters 2 lockable brakes
- › Lifting range 1310 mm
- › Lifting capacity up to 440 pounds lbs. or 200 kg
- › Emergency on/off button
- › Manual and electric emergency lowering
- › Removable 24 V / 2,9 Ah Battery-Pack incl. wall charging device
- › Comprehensive range of slings to meet your individual requirements
- › Color: Pearl Gold / Anthracite
- › The LEXA patient lift is also available with a scale option, to record patient weights



A multi-talent for careful transport of patients

The Original LEXA bariatric battery-powered lift is an absolute multi-talent among lifts. It is a very comfortable belt lifter, which, due to its swivel spreader bar, has many capabilities and applications. The lift offers motorized patient positioner and leg spreading and is available with either a wide variety of lifting slings.

For instance, it can easily take up patients sitting, lying, or from the floor and move them easily and safely to a sitting or lying position so they can be transported gently and comfortably to the toilet, to the wheelchair, or to the bed. The sway-free transport also conveys a feeling of security.

Everything can be done with a user-friendly electronic hand held control that allows the caregiver to properly focus on patients. Unique Safety Features make the LEXA Patient Lift stand apart from other.



Base length	1220 mm
Base width closed/open	668/1050 mm
Base height	117 mm
Overall height min./max.	1750 mm / 2240 mm
Lifting range	1310 mm
Patient positioner	4-point clip connector
Lifting capacity	up to 440 pounds (lbs.)
Battery capacity	24 V / 2.9 amps

Barrier Free Lifts, Inc. | 1620 SW 17th St | Ocala, FL 34471

Phone: 1-800-582-8732 | Fax: 1-866-378-3318 | us-office@barrierfreelifts.com | www.barrierfreelifts.com

**BARRIER
FREE LIFTS**