

Day Chair range

Air comfort



The Ultimate in
PRESSURE CARE



K Care

kcare.com.au

Deluxe Chair V2 - **Small** [AC59114](#)



Deluxe Chair V2 - **Medium** [AC59116](#)



Deluxe Chair V2 - **Standard** [AC59111](#)



Deluxe Chair V2 - **Maxi** [AC59115](#)



Day Chair features

The Air Comfort Deluxe V2 is the market leader in pressure care seating. With unparalleled pressure redistribution and adjustable positioning, the Air Comfort Deluxe V2 offers ease of handling for carers and maximum comfort for patients.

- The Patented *Airgonomic Seating System* uses four air cells to immerse and cushion the patient, providing comfort and pressure redistribution.
- Adjustable back recline, seat tilt and raising leg rest for patient comfort and positioning.
- Drop away /adjustable arm rest's and back wings for improved access.
- 125mm (150mm Maxi) Twin wheel Tente castors (Three total lock and one directional) for safe handling and manoeuvrability.
- Carrflex Extreme multi-stretch material with *Sanitized* treatment and enhanced durability.
- Access Vehicles Approved (Australian Standards) for transport in Access Vehicles.
- Replaceable parts and upholstery cover for easy cleaning and maintenance.

Day Chair accessories



Meal tray - Small [AO59036](#)
 - Large [AO59035](#)
 - Maxi [AO59007](#)



Incontinence cover - Std [AO59013](#)



Seat belt buckle [AO59014](#)



Foot rest - Std [AO59055](#)
 - Maxi [AO59053](#)

Dimensions					
	Overall	Internal seat	Recommended patient height	Product weight	Safe working load
Small	L: 1000mm W: 670mm H: 1130mm	W: 430mm D: 400mm H: 550mm (from floor)	Up to 1420mm (4ft 8")	40kg (approx.)	180kg
Medium	L: 1100mm W: 670mm H: 1130mm	W: 430mm D: 470mm H: 550mm (from floor)	1420 - 1620mm (4ft 8" - 5ft 4")	40kg (approx.)	180kg
Standard	L: 1130mm W: 680mm H: 1260mm	W: 470mm D: 550mm H: 570mm (from floor)	1620 - 1880mm (5ft 4" - 6ft 2")	45kg (approx.)	180kg
Maxi	L: 1360mm W: 950mm H: 1385mm	W: 720mm D: 550mm H: 640mm (from floor)	1620 - 1880mm (5ft 4" - 6ft 2")	100kg (approx.)	300kg

Please note: All specifications shown are subject to change without notice

Day Chair options

Electric - Medium [AC59316](#)
-Standard [AC59112](#)

Full Electric / Battery operation is available with Linak systems to operate the back recline, seat tilt and raising legrest functions.



Alternating Pressure Care Overlay - Standard [AC59441](#)

A battery powered alternating overlay can be fitted to the Deluxe V2 for full active pressure care.

