

Technical Data


Electrical data, dimensions and weights

	DIANA Comfort	DIANA Standard
Net weight (without Battery)	33,5 kg	32,5 kg
Base length	1187 mm	1160 mm
Base width (outside)	683 - 1043 mm	670 - 1029 mm
Base width (inside)	527 - 887 mm	540 - 899 mm
Base height	137 mm	107 mm
Chassis clearance	75 mm	45 mm
Overall height	1446 - 2013 mm	1416 - 1983 mm
Lifting range	989 mm	989 mm
Power supply	24 V	24 V
Battery capacity:	24 V / 2,9 Ah	24 V / 2,9 Ah
Duty cycle	10% (1:9)	10% (1:9)
Degree of protection	IPX5	IPX5
Noise Level		
Idle	< 48 dB	< 47 dB
Max.	- 48 dB	- 47 dB
Battery charger:	Wall mounted charger 230 V / 400 mA IPX4	Charging cable 230 V / 260 mA IPX4
Lifting capacity	200 kg / 440 lbs.	160 kg / 350 lbs.

Technical modifications reserved.

Technical Data

Dimensions DIANA Comfort



Technical Data

Dimensions DIANA Standard

