

Technical Data

Electrical data, dimensions and weights

	MONA
Base length:	815 mm
Base width:	990 mm
Overall height:	1350 - 2050 mm
Floor clearance:	70 mm
Base height:	155 mm
Lifting range:	700 mm
Power supply:	36 V
Battery capacity:	36 V / 2.3 Ah
Battery charger:	Horcher Series 100 Typ IP20 T 35/E 120 V/ 36 V/0,4 A
Fuse type for charger	T 0,63 A in charger Type IP 20 T 35/E, Type G 136-0,4 , Horcher
Current for loading station:	5 A
Noise level:	
Idle	< 49 dB
Max.	- 49 dB
Net weight	70 kg / 154 lbs.
Lifting capacity	160 kg / 350 lbs.

Technical modifications reserved.

Technical Data

Dimensions

