



Specifications

Electrical data, dimensions and weights

	Raisa-Pro
Net weight (without Battery)	80 kg
Lifting capacity	275 kg / 600 lbs.
Lifting range	645 mm
Base length	1239 mm
Base width min. / max.	670 / 1333 mm
Base height	132 mm
Overall height / Mast	1427 mm
Overall height / Lifting arm	1427 / 1933 mm
Chassis clearance	56 mm
Caster diameter	100 mm front, 125mm rear
Turning radius	550 mm
Sling suspension	2-point-clips
Noise level	
Idle	< 44 dB
Max.	- 44 dB
<u>Electric</u>	
Power supply	24 V
Battery capacity:	24 V / 2,9 Ah
Duty cycle	10% (1:6)
Fuse	T 6 A
Degree of protection	IP54
<u>Battery charger</u>	
Series	Typ G124-1 T 35/E
Prim. / Sek. / I max.	120 V ≈ / 24 V = / 1 A
Sicherungstyp	T 0,63 A
Schutzart	IP20
Schutzklasse	II

Technical modifications reserved.

Specifications



Specifications

Dimensions RAISA-Pro

