



Specifications

Electrical data, dimensions and weights

	LEXA-Pro
Base length	1239 mm
Inner base length	773 mm
Base width closed/open	1333 mm / 668 mm
Inner base width closed/open	1203 mm / 596 mm
Floor clearance	33 mm
Base height	102 mm
Overall height min.	1757 mm
Overall height max.	2217 mm
Lifting range	1303 mm
Height of CSP (max./min.)	1678 / 375 mm
Range of CSP	702 mm
Min. Turning Radius	ca. 665 mm
Net weight	ca. 70 kg
Lifting capacity	275 kg
Noise Level	
Idle	< 44 dB
Max.	- 44 dB
<u>Electricity</u>	
Power supply:	24 V
Max. Power Consumption	120 W
Battery capacity:	24 V / 2,9 Ah
Duty cycle	10% (1:6)
Fuse	T 6 A
Application part	 Typ B
Degree of protection	IP54
<u>Battery charger</u>	
Model	Type G124-1 T 35/E
Prim. / Sec. / I max.	230 V ≈ / 24 V = / 1 A
Fuse	T 0,315 A
Degree of protection	IP20
Protection class	 II

Technical modifications reserved.