



Healthcare Equipment

K Care Healthcare Equipment

100 Mulgool Rd
WA 6090, Australia.
Telephone +61 8 9248 4444
Facsimile +61 8 9248 4190
www.kcare.com.au

**APPROVED CLEANING
PRODUCTS FOR MATTRESS
COVERS, UPHOLSTERY
TOPS AND SHOWER
TROLLEY TUB FABRICS:**

CLEANING PRODUCT:	MAXIMUM CONCENTRATION:
Clinell Hand & Surface Sanitizer Disinfectant Clinell Universal Sanitising Wipe Clinell Detergent Wipes Clinell Sporidical Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
PAL Detergent Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
Mediwipes Disinfection (Sulco Limited, New Zealand)	Use according to the manufacturers specifications.
PAL Disinfectant Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
Surface Sanitising Wipes Neutral Detergent Wipes (Reynard Health Supplies, Australia, New Zealand & Europe)	Use according to the manufacturers specifications.
PureGreen 24 (Pure Green 24, America)	Use according to the manufacturers specifications.
Spring Clean (Eco Lab, Australia)	5% Spring Clean (1 part Spring Clean to 20 parts water 50ml per litre)
Ecopine (Worldwide)	5% Ecopine (1 part Ecopine to 20 parts water 50ml per litre)
Actimax (Worldwide)	Use according to the manufacturers specifications.
Viraclean (Whiteley Medical, Australia)	Use according to the manufacturers specifications.
Virkon (Dupont, Worldwide)	5% Virkon (1 part Virkon to 20 parts water 50ml per litre)
Trigene (Bio Strategy, Australia & New Zealand)	Use according to the manufacturers specifications.