

Portable Ceiling Lift UPL / PC-2



The advantages at a glance

Portable ceiling lift for more freedom of motion

- › Can be used in wet and dry rooms
- › Light and compact with only 8,5 kg / 19 lbs.
- › Easy, smooth and safe transport
- › Easy hook up and unhook
- › 2 functions by hand control
- › Electric emergency down function
- › IBS - Intelligent Battery System
- › 36 Volt / 3.0 Ah Battery Pack
- › Lifting speed 6-8 feet per minute
- › Can be used in different track systems
- › Comprehensive range of slings to meet your individual requirements
- › Lifting capacity 600 lbs. / 275 kg
- › Available with clips or loops suspension

We all want to be able to move about safely, in comfort and with dignity and the Horcher PC-2 has proven itself to be an invaluable aid to thousands of people with restricted mobility. Capable of being installed in virtually unlimited configurations, from simple track systems to room covering transversals, our technical engineers will be happy to design track layouts to meet your individual requirements. In whatever location, whether at home, in the workplace, near the swimming pool, or almost anywhere you desire, PC-2 can meet your needs!

The PC-2 can easily be taken off one track, carried to another, and reattached, thus avoiding the need to purchase more than one lift for several different locations. Furthermore, installing small pieces of track where you actually need it saves on more costly continuous track systems. Available with the most comprehensive product line of slings and models with lifting capacities of up to 600 lbs, the Horcher PC-2 just may be the lift solution that you have been looking for.



So reliable and easy to use!

Suitable slings for the UPL / PC-2 lift system

